

Biocraft's Electronic Test and Measurement Instruments

Economical but Tested, Guaranteed Best Equipments

CEI-1 **'Craft's Digital Ph Meter** Complete with combined Electrode, Buffer solution stand and working manual. Range 0-14 pH with auto temperature compensation 3.5 digital led display

Rs. 4,800

CEI-2 **'Craft's Pocket pen type PH meter** no need of electrode. Instant Reading Imported

Rs. 1,350

CEI-3 **Craft's MicroProcessor Bases pH Meter**—with Alpha-numeric LCD display showing pH, mV, Temp. having auto temperature compensation, auto buffer recognition with all accessories and manual

Rs. 15,200

CEI-6 **'Craft' Digital Photoelectric Calorimeter** with 8 filters with facility for automatic zero setting range 400 to 700 nm.

Rs. 7,200

CEI-7 **Craft Lux meter** for photometric measurements. Complete with standard accesories

Digital Lux Meter—0—50,000 Lux IMPORTED

Rs. 3,590

CEI-8 **Craft LUX meter** 0-2 lacs LUX

Rs. 7,900

CEI-9 **Craft Turbiditimeter/Nephelometer**

A Digital Range upto 1,000 NTU/JTU

Rs. 9,550

CEI-10 **Craft Digital T.D.S. meter** (Total DISSOLVED solids) complete with accessories

Rs. 6,200

CEI-11 **Craft Digital Potentiometer** with reference & Platinum electrode. Range 0—2000 mV

Rs. 6,100

CEI-14 **'Craft' Digital Conductivity meter** with manual

Rs. 5,700

CEI-15 **Dissolved oxygen Meter** (Digital) with accessories

Rs. 7,300

CEI-16 **Glucometer** Electronic digital will accessories, Manual and Testing strips


Rs. 3,000

CEI-17 **Digital Spectrophotometer (Double Display)**—Range 340-960 nm with facility for automatic concentration, % Transmission and adsorbance measurements. 3-1/2 digit display, separate 3 digit LED display for Wavelength resolution. Equipped with a set of 2 Cuvettes.

Rs. 30,800

CEI-17A **Microprocessor Spectrophotometer Visible Single Beam**—Range 320-1020 nm facility for Auto concentration, % Transmission, Absorbance and K-Factor measurement, It has 3-1/2 Bright LCD display for data. Its sample holder accommodates glass cuvettes of 10 to 50 mm path length. Facility has been provided for unlimited sample data storage and interface for any printer through RS232C. A set of 4 glass cuvettes of 10 mm path length is provided. Including software.

Rs. 42,800



Biocraft's Electronic Test and Measurement Instruments

CEI-17B **Single Beam Microprocessor UV-Vis Spectrophotometer (Advance Model)**—Range 195 to 1000 nm with facility for automatic concentration, % T, Absorbance & K Factor measurement. Automatic Wavelength calibration and selection. It has sample holder with 4 position spectro control for 10 mm path length cuvettes. Facility has been provided for interface with any printer USB. Set of 2 quartz and 4 glass cuvettes is provided. Highly stable and accurate. Software included. (Imported)



CEI-17B

Rs. 1,55,000

CEI-17D **Double Beam Microprocessor UV-Vis Spectro-photometer**—Range 190-1100 nm, fully automatic operations, automatic 8 cell change, bandwidth 1 nm, automatic wavelength calibration & wave length scan, automatic peak picking & spectrum display. Large screen with graphic display. Spectrum and data can be printed out by printer and set to computer via USB port. Set of 4 quartz and 16 glass cuvettes is provided, software included. (Imported)




CEI-17D

Rs. 3,71,000

CEI-18 **Auto Karl Fischer Titrimeter** with Separate control and Stirrer unit, automatic burette, dual platinum electrode & accessories. **Rs. 18,800**

CEI-19 **Craft's Digital Salinity Meter** with Salinity Cell & Accessories. **Rs. 7,400**


CEI-20 **Digital Flame photometer (Single Channel)** with filter (Na & K) 3-1/2 digit LED display (ppm) along with compressor and other accessories, auto ignition, Ca & Li filter can be provided on request at extra cost.



CEI-18

Rs. 32,200

CEI-21 **Microprocessor Flame Photometer**—Four Element measurement in a single aspiration (Na, K, Li and Ca) Calibration curve fit routines first order and second order, Data displayed in concentration units, i.e., ppm & meq, automatic filter selection; auto calibration feature. Printer attachment, PC connecting compatibility, Sample data storage more than 500 results. Alphanumerical key board for easy data entry. Graphical LCD display thus easy to read, Curve calibration using upto 5 standard, Multiple calibration curve saving facility, stainless steel burner. Auto gas cut off in case of power failure provided with 2 filters and compressor. **Rs. 65,200**



CEI-20

CEI-23 **Digital Haemoglobinometer**—Microcontrolled based square cuvette STD and non STD model sample size 1 ml **Rs. 9,950**



CEI-21

ELECTRONIC POCKET TESTERS IMPORTED

CEI-35 Pocket PH meter 0-14 pH/0.2 pH **Rs. 1,350**

CEI-36 T.D.S. Tester 10-1990 PPM Digital **Rs. 1,400**

CEI-37 T.D.S. Tester 100-10,000 PPM Digital **Rs. 1,400**

CEI-38 Conductivity Tester Digital 10-1990 US **Rs. 1,700**


CEI-39 Conductivity Tester Digital 100-19,900 uS **Rs. 1,700**

CEI-40 Salinity Tester Digital 0.00-1.00 P NaCl **Rs. 8,400**

CEI-43 Soil pH and Moisture Meter Take Mura **Rs. 3,850**

CEI-44 Digital Sound Level Meter Imported Range 30 to 130 dB (3 Ranges) **Rs. 8,400**

CEI-45 Grain Moisture meter LSKC-60 (Paddy, Rice, Corn, Wheat) **Rs. 8,200**



POCKET TESTER