


## Biocraft's Microscope Accessories

**BMA-1 Double demonstration eye-piece with imported prism**—Prismatic convenient for simultaneous viewing enabling the teachers to demonstrate the students. Complete in velvet lined case. **Rs. 1,100**



**BMA-1**

**BMA-4 Camera lucida (Prism type)**—Imported prism but assembled in India to fit shoulder of the microscope. Complete in a velvet case. **Rs. 900**

**BMA-6 Camera lucida (Mirror type)**—Sturdy and reliable equipment fits on the eye tube of all microscopes and has large mirror with adjustable inclination. Complete in a velvet lined case. **Rs. 1,380**

**BMA-8 Microscope lamp (Substage)**—(Pin type model)—(Special quality)—Can be fixed under the stage of microscope with LED white light. Battery operated. **Rs. 440**

**BMA-10 Craft Laboratory Table Lamps**—sturdy and stream lined body Angle Poised.

(1) One Arms, Extension Approx. 14"


**Rs. 415**

(2) Two Arms Extension Approx. 24"

**Rs. 1,180**



**BMA-4**



**BMA-7**



**BMA-8**



**BMA-10**

**BMA-11 Optic Set**—Consist of two Huygenian Eyepieces, 10x, 15x and two achromatic objective 10x, 45x. (German lenses assembled)

**Rs. 1,250**

**BMA-12 Optic Set**—Same as above Imported make

**Rs. 3,800**

**BMA-13 Oil Immersion Lens 100x**—Japan/Meopta/German/Imported

**Rs. 2,700**

**BMA-14 Oil Immersion Lens 100x**—German/Japan Lenses assembled in India spring loaded

**Rs. 1,450**

**BMA-15 Oil for Oil Immersion**—**Repacked** Imported 100 ml.

**Rs. 185** Per bottle

**BMA-16 Eye Piece (Imported)**—5x or 10x or 15x

**Rs. 250** each

**BMA-17 Eye Piece (Indian)**—Assembled with imported Lenses 10x, 15x.

Each **Rs. 175**

**BMA-18 Eye Piece wide field (WF)** 10x or 15x mm coated

**Rs. 295**

**BMA-20 Eye Piece for dissection microscope** 10x or 20x

**Rs. 55**

**BMA-21 Spare Objective**— 2x 3x or 5x 10x 45x

Imported

—

**Rs. 750**

**Rs. 1,180**

Assembled (Indian)

**Rs. 250**

**Rs. 320**

**Rs. 490**

**BMA-22 Reflector Mirror** for Microscope Complete with Steel Frame




**Rs. 95**



**BMA-21**


**Biocraft's Microscope Accessories**

- BMA-29 Microscope **Stage Clips** **Rs. 60** pair
- BMA-30 **Pointer Eye Piece Imported** 15x, 10x, 5x **Rs. 370** each
- BMA-31 **Pointer Eye Piece Indian** 10x, 15x, 5x **Rs. 295** each
- BMA-33 **Micrometer Eye Piece Imported**—10 mm divided in 100 div. **Rs. 425**
- BMA-35 **Graph Eye Piece Imported**—One square cm. is divided in 400 Squares **Rs. 410**
- BMA-37 **Stage micrometer objective imported** 1mm divided into 100 (0.01 mm) **Rs. 540**
- BMA-38 **Stage Micrometer objective imported** 10 mm divided into 100 (0.1 mm) **Rs. 310**
- BMA-39 **Stage Micrometer (Indian)** **Rs. 250**
- BMA-40 **Ocular micrometer disc**—  
 A—10 mm divided into 100 division—**Imported Rs. 210**  
 —do— **Indian Rs. 150**
- BMA-41 **Ocular graph disc**—Imported 0.5 × 0.5 × 10 mm **Rs. 260**
- BMA-43 **Ocular cross line disc** **Rs. 250**
- BMA-44 **Lens Cleaning paper**—100 sheets, per book **Rs. 70** per book
- BMA-45 **Lens cleaning solution** per bottle **Rs. 60**
- BMA-46 **Mechanical Stage**—Attachable to any standard microscope. All Brass. **Rs. 1,400**
- BMA-50 Plastic (thick) covers for Microscopes **Rs. 48**
- BMA-51 Plastic covers for Binocular Microscope **Rs. 80**
- BMA-52 Plastic cover for Microtome **Rs. 195**



**Biocraft's Mini Cleaning and Repair Kit**

- MCR-1 **Mini Cleaning and Repair Kit**—It can be set right than and there with 'Mini Kit' containing. (1) Long nose pliers; (2) Screw driver (set of 6); (3) Small fine forceps 4" or 5" SS; (4) File triangular; (5) A bottle for oiling; (6) Lens cleaning paper; (7) Lens cleaning solution; (8) Sand paper, Flannel, Napkin, Brush. Complete kit in a Plastic box **Rs. 710**
- MCR-2 **Delux-Mini Cleaning and Repair Kit**—same as BMPK-1 but with more superior quality tools made by Taparia or equivalent, known for their everlasting service. It is supplied with soldering iron. **Rs. 1100**
- MCR-3 **LENS CLEANING KIT** : Ideal for cleaning lens of Microscopes, camera & other optical Instruments in box **Rs. 450**
- MCR-4 **TOOLS KIT** (Imported) having 100 pieces for all type of repairwork **Rs. 3,100**



## Biocraft's magnifiers

BMF-1 **Craft's Magnifier folding**—Imported lenses, mounted in India in plastic frame.

Single (4x) **Rs. 48**  
 Double (7x) **Rs. 52**  
 Triple (10x) **Rs. 78**

BMF-2 **Craft's Magnifier folding**—Gowland England type in metallic frame 10x suitable for field collection. With fine thread for hanging **Rs. 60**

BMF-4 **Imported Hand magnifiers** with classic design and large viewing Area with handle

3" — **Rs. 180**      4" — **Rs. 220** each

BMF-5 **Craft's Magnifier with handles**—

2" — **Rs. 49**      3" — **Rs. 110**      4" — **Rs. 160**

BMF-6 **Craft's Magnifier torch flash magnifier**—With magnifier lens attached 13x or 7x **Rs. 270**

BMF-7 **Craft's Megascop**—Fitted with 6" magnifier and round florescent tube and transformer. Very useful for lab use. **Rs. 4,450**

BMF-8 **Craft's Magnifier on stand**—With heavy base nickle plated rod on which the magnifier horizontally moves up and down.

3" — **Rs. 590**      4" — **Rs. 740**

BMF-9 **Craft's Dissecting magnifier**— (Horse shoe heavy base)—Lens adjustable to any direction.

3" — **Rs. 750**      4" — **Rs. 980**

BMF-10 **Craft's Crane Magnifier**—Long flexible arms enables the lens to be adjusted to any distance and direction

3" — **Rs. 750**      4" — **Rs. 980**



BMF-11 **Wide Field Magnifier**— **Rs. 160**




BMF-1



BMF-2




BMF-6




BMF-7



BMF-11



BMF-10



BMF-9

## GUARANTEE

ALL BIOCRAFT'S Instruments are covered with ONE YEARS GUARANTEE' against any manufacturing defect please call us if you face any trouble with our instruments.