







Biocraft's Microtome Accessories

- BMTA-1 **Craft's Knife for rotary micro-tome with handle**—Manufactured from fine quality carbon steel with back and handle in box.
 120 mm 140 mm 160 mm 180 mm
Rs. 1500 Rs. 1700 Rs. 2000 Rs. 2,450
- BMTA-2 **Craft's Imported Microtome Razor.** Stainless steel 120 mm
Rs. 4,600 each
- BMTA-3 **Craft's Knife for rocking micro-tome with handle.**
Rs. 900
- BMTA-5 **Craft's Section Blade**—(Imported for blade holder)
Rs. 60 each
- BMTA-6 **Craft's 'L' mould for casting paraffin block**— Made of brass.
 Small size 37 × 25 × 15 mm **Rs. 250 pair**
- BMTA-7 -do- Large size 50 × 25 × 15 mm **Rs. 290 pair**
- BMTA-8 **Craft's Imported Steel structure.** For making many blocks at a time
Rs. 750 each
- BMTA-9 **Object disc**—(Block holder)—For Rotary microtome, set of 3
Rs. 300
- BMTA-12 **Craft's Embedding Rings**—Disposable plastic embedding moulds 25×25×12 mm
Rs. 48 doz.
- BMTA-13 **Craft's tissue capsule stainless steel** with perforation and sliding cover very helpful during the process of dehydration. 16×16×12, 22×22×12, 32×25×12 mm
Rs. 100
- BMTA-13 A. **Craft's Embedding Cassette** (Disposable) plastic
Rs. 10 each
- BMTA-16 **Craft's Sharpening Hones**—Made of specially selected uniform ingredient to provide convenient, rapid and satisfactory sharpening.
 6" × 2" — **190** 8" × 2" — **285** 12" × 2.5" — **1,550**
- BMTA-19 **Diamond Point Glass Marking Pencil** **Rs. 95**
- BMTA-20 **Soft (Coloured) Glass Marking Pencil** **Rs. 18 each**
- BMTA-22 **Slide Ringing Table**—A turn table for slide ringing Attractive finish, smooth running with ball bearing
Rs. 1,520
- BMTA-24 **Slide Ringing Solution**—Black or transparent per 50 ml.
Rs. 110
- BMTA-26 **Craft's Staining Assembly**—with hanger
 Capacity Aluminium
 25 slides **Rs. 175**
 50 slides **Rs. 270**
- BMTA-27 **Staining Square Bar with stand**—35 × 65 cm
 Aluminium **Rs. 198**
 Stainless steel **Rs. 265**




Biocraft's Microtome Accessories


- BMTA-29 Staining Slide Tray Aluminium** (Double purpose)—Useful for staining slides in flat position or drying freshly prepared slides in vertical position. **Rs. 130**
- BMTA-36 Imported Disposable Microtome Blades** pack of 50 pc in dispenser. It dispenses one blade at a time. The blade is made of stainless steel size 80 mm length. Low profile or high profile. **Rs. 8,900** per packet
- BMTA-37 Craft's Disposable Blade Holder**—Suitable for holding low and high profile blades universal design and can be filled on all type of Rotary Microtomes. **Rs. 6,500**
- BMTA-38 Craft's Tissue Embedding Moulds**—Made from stainless steel. Eliminates need for trimming of blocks.
- | Size in mm | Price |
|--------------|---------------------|
| 16 × 16 × 12 | Rs. 115 each |
| 22 × 22 × 12 | Rs. 115 each |
| 32 × 25 × 12 | Rs. 115 each |
- BMTA-40 Craft's Paraffin Dispenser**—Very useful for Histo pathology and Microtome labs for making wax embedded mould or Blocks A-Tap with Temperature Control is provided and the wax flows freely thermostatically controlled unit. Capacity—5 Lt. **Rs. 16,900**



BMTA-29



BMTA-36



BMTA-37



BMTA-38



BMTA-40

