

Biocraft's Animal Physiological Apparatus and Accessories

BPZA-1 Craft's Electric Recording Drum (Kymograph)—It is made on Polmar pattern. Double electric contract arm and interchangeable contract blocks are provided as standard fixture. The gear train comprise of two mill cut gears.


Speed range from 0.12 to 640 m.m. per second and is accurate to $\pm 1\%$. The change in the speed from one to the other can be effected through a lever adjustable, provided with electric long duty motor and the speed of it is controlled by mechanical governor control device to ultimate rate of 1 : 2 *Electric Drum is supplied with Chromium plated steel spindle and 6" x 6" cylinder with screw lift.*

Rs. 9,200

BPZA-3 INCO Electric Recording Drum—E/8 Heavy duty motor with 8 Ratio gear box complete with stainless steel spindle.

Rs. 16,800

BPZA-6 Craft's Physiology Kit No. 1 smoke writing. The kit includes—



Sl.No.	Quantity	Description	Rs. (each)
1	1	Kymograph—BPZA—1	9,200
2	25 sheets	Kymograph—paper 6 x 20"	120
3	1	Induction coil with variable interrupters	1,950
4	1	Adjustable stand with 15" up right.	820
5	1	Muscles lever (Bodies universal liver)	210
6	1	Weight and Hook Set 11 gm.	210
7	1	TAMBOUR—Brass-bright CP.	340
8	1	Tuning fork steel N 100	270
9	1	Tuning fork starter.	40
10	1	Femur clamp of brass bright CP.	580
11	2	Double clamps (X block brass)	400
12	1	Electric time maker.	810
13	1	Heart lever brass.	310
14	1	Frog board with clamp and plastic bands.	380
15	1	Stylus with pointer Aluminium	60
16	1	Simple key with ebonite base.	290
17	1	Electrode with copper wire	280
18	1	Smoking burner gas.	1,250
19	1	Muscles chamber with muscles lever and electrodes.	1,050
20	1	Heart chamber with heart lever and electrodes.	1,050
Price of full kit.			Rs. 19,670

BPZA-12 STUDENT ORGAN BATH—Made of Acrylic thick sheet complete stand, levers, oxygen tube and tubular heater and stirrer

Rs. 4,100

BPZA-12 A. Student Organ Bath—Double Unit

Rs. 8,800

BPZA-13 Student Electronic Stimulator

Rs. 2,850

BPZA-16 Warberg's animal respirometer—A glass apparatus to experiment on respiration in animals.

Rs. 700 each