

ELECTRONIC WATER AND SOIL ANALYSIS KITS

- CEI-60 **Craft's Portable water/Soil Analysis kit**—Digital 3 Parameters measures pH, conductivity and ORP (Oxidation Reduction potential) **Rs. 8,900**
- CEI-62 **Craft's Portable water/Soil Analysis kit**—Digital 4 Parameters measures PH, ORP, Conductivity and Temperature **Rs. 11,900**
- CEI-64 **Craft's Portable water/soil kit**—Digital 5 Parameters—PH, ORP, conductivity, Temperatura and DO, (Dissolve Oxygen) **Rs. 14,500**
- CEI-66 **Craft's Portable water/soil kit**—Digital 6 Paremeters—PH, ORP, conductivity, Temp, Dissolve Oxygen and Tota Dissolved Solids. **Rs. 15,400**
- CEI-68 **Craft portable water/soil kit**—Digited 7 Parameters (+ Salinity) **Rs. 16,850**
- CEI-70 **Craft's portable water/soil kit** with 8 Parameters : pH MV conductivity, TDS Salinity, Turbidity, D.O. and Temp. **Rs. 22,000**



Biocraft's Soil Testing Equipment

- SOL-1 **Apparatus for Soil Determination**—consisting of a core cutter 100 mm dia dolly and a rammer with steel rod. **Rs. 4,300**
- SOL-2 **Cone Penetrometer**—for determining the liquid limit of soils **Rs. 3,600**
- SOL-4 **Sampling Auger outfit**—posthole type complete with one meter rod, and handle etc.
 38 cm **Rs. 2,750** 50 cm **Rs. 3,200** 75 cm **Rs. 3,480**
- SOL-6 **Soil moisture sample box** (S.S.) Pack of 10 boxes **Rs. 520 per set**
- SOL-8 **Soil Keen's box** for soil analysis **Rs. 3,100**
- SOL-9 **Pocket Penetro Meter** **Rs. 4,200**

Biocraft's Accessories for Electronic Instruments

1. Combined Electrode **Rs. 780** each
2. Glass Electrode **Rs. 800** each
3. Platinum Electrode **Rs. 1,250** each
4. Gelfilled Compound Electrode **Rs. 1,300** each
5. Silver electrode **Rs. 1,100** each
6. Reference Electrode **Rs. 700** each
7. Square Cuvettes—Glass—950, Quartz **Rs. 2,100** each
8. Conductivity cell (0.5) **Rs. 2,450** each
9. Temperature Probe **Rs. 800** each
10. Selenium Photo Cell **Rs. 1,450** each
11. Standard Glass Filler **Rs. 300** each
12. Cuvettes (10 ml) glass—set of two **Rs. 950** each
13. Cuvettes (10 ml) Quartz **Rs. 2,650** each

