

Biocraft's Entomological Accessories

BWEA-1 Craft's Insect Boxes—(Entomological Boxes) used for the storage of pinning insects collection. The bottom of the case is lined with P.U. sheets. Boxes are proved with groves for keeping Insecticides

18"×12"×3" 15"×12"×3" 12"×9"×2½"
Rs. 720 **Rs. 550** **Rs. 495**



BWEA -1

BWEA-2 Craft's Insect Pinning Case (Deluxe Type)—Same BEWA-1 out with fine-wood and finish & lining is of imported EVA sheet.

18"×12"×3" 15"×12"×3" 12"×9"×2½"
Rs. 1,390 **Rs. 1,150** **Rs. 710**



BWEA -2


BWEA-3 Craft's Insect or dry plant display case with glass top—Double frame glass top. The bottom of case is lined with P.U. Sheet.

18"×12"×3" 15"×12"×3" 12"×9"×3"
Rs. 790 **Rs. 625** **Rs. 570**

INSECT TRAPS AND INSECT CABINETS

BWEA-5 Craft's Insect Night Berlese Funnel—Large metallic funnel, circular with light shade electric fittings. Very useful for night collection and soll arthropods. Complete with stand.

Rs. 1,400



BWEA -5

BWEA-7 Craft's Insect Suction Trap—For collecting small insects with long corning glass tube fitted with Latex tube and Brass gauge to protect the entry of insects in the mouth.

Rs. 210



BWEA -7

BWEA-8 Craft's Insect Pooter—A suction trap specially designed for small Diptera and other insects. Complete with glass vial and rubber tubes.

Rs. 210

BWEA-9 U.V. Insect Trap—Equiped with ultraviolet imported tube very useful for night flying insect collection.


Rs. 1,680



BWEA -9
U. V. Insect trap

BWEA-10 Craft's Metallic Insect Light Trap—Madea mild steel size 40×40×25 cm with 5 cm dia metal funnel coming into the box for convenient collection of trapped insects.

Rs. 7,950



BWEA -10

BWEA-11 Craft's Insect showcase cabinet—Made of board frame work covered with thick ply. Single door which dust and pest proof. Showcase are finished in stain and moisture proof nitro-cellulose clear lackuer sliding on wooden rail. Lined with P.U. Sheet

(A) **Cabinet with 6 showcases** of Size—
12" × 10" × 2½"

Rs. 7,450

(B) **Cabinet with 12 showcases** of Size—
12" × 10" × 2½"


Rs. 10,500

(C) **Cabinet with 6 showcases** of Size—
18" × 12" × 3"

Rs. 9,100

(D) **Cabinet with 12 showcases** of Size—
18" × 12" × 3"

Rs. 12,400



BWEA -11

Biocraft's Entomological Accessories

BWEA-11-A. **Craft's Insect Showcase Cabinets**—Super delux quality with superior wood fitted with 'EVA' sheet for pinning the insects.

(A) **Cabinet with 6 showcases** of Size—12" × 10" × 2½" **Rs. 8,850**

(B) **Cabinet with 12 showcases** of Size—12" × 10" × 2½" **Rs. 11,900**

(C) **Cabinet with 6 showcases** of Size—18" × 12" × 3" **Rs. 13,950**

(D) **Cabinet with 12 showcases** of Size—18" × 12" × 3" **Rs. 16,100**

BWEA-12 **Craft's Insect collection bag**—Made of thick canvas cloth with adjustable strap to carry 6 collection bottles. **Rs. 480**

BWEA-13 **Craft's Insect pinning forceps**—S. S. curved and blunt. **Rs. 58**

BWEA-14 **Craft's Pinning block**—Made of Alluminium block to fix the insect at the desire level of insect pin. **Rs. 250**


BWEA-15 **Craft's Entomological dissection box**—Set of 8 stainless steel fine instruments including teasing needle made of insect pin. **Rs. 390**

BWEA-16 **Craft's Entomological Collection Kit**—The kit include insect net, pinning forcep, plastic Tubes—6 Nos., Insect preservative Killing Jar, P.U. Sheets, (6 Nos.), Paper envelops. Full kit in a carrying Bag. **Rs. 2,150**

BWEA-17 **Craft's Integrated Pest Management Kit (IPM)**—Kit of 14 items with Insect Collection Net **Rs. 1,980**



BWEA-12



BWEA-18

CRAFT'S ENTOMOLOGY PINS

BWEA-18 **Craft's Entomology pins**—English or Austria made from fine tempered steel with needle points and secure heads black or white Minutent, 00, 0, 2, 3, 5, 12, 16, 20 Per pkt. of 100 Nos. **Rs. 450**



BWEA-20

INSECT KILLING BOTTLE INSECTICIDES AND INSECT CAGES

BWEA-19 **Craft's Insecticide powder**—Very strong insecticide for keeping in insect boxes and dry collections to protect them from insects. per 500 gms. **Rs. 172**

BWEa-20 **Craft's insect killing jar**—Screwed bakelite cap, wide mouth jar charged with ethyl acetate. **Rs. 148**

BWEA-21 **Craft's Insect preserving solution**—A standard preservative. 500 ml. **Rs. 195**

BWEA-22 **Paradichloro Benzene**—500 gm **Rs. 380**

BWEA-23 **Craft's Insect spreading board (adjustable)**—Size 14"×5" Groove adjustable with metallic fittings with cork sheet. **Rs. 295**

BWEA-24 **Craft's Fixed spreading board**—Non-adjustable with cork superior **Rs. 240**


BWEA-26 **Craft's Insect Rearing cage**—Made of aluminium framework with wire mesh on three sides. Provided with sliding or hinged door. Very useful for field collection as well as laboratory use.

Size 6"×6"×6"


Price **Rs. 350**

8"×8"×8"

Rs. 370



BWEA-23



BWEA-26