

Biocraft's Glass Apparatus Made of Borosilicate Glass

DA-1 **Micro kjeldal distillation apparatus**—Pregl's for determination of Nitrogen automatic discharge according to 'Parnes and Wanger' made of Corning Glass with silver tube, clamp stand and burner. **Rs. 11,900**

DA-2 **Distillation Assemblies (Quick fit type)**

I. (a) R.B. Flask—2000 ml, (b) Splash Head, (c) Double Surface condenser, (d) Receiver adopter. Complete unit **Rs. 2,700**

II. (a) R.B. Flask—1000 ml, (b) Recovery bend, (c) Thorpe Condensor, (d) Receiver adopter **Rs. 3,000**

DA-3 **Soxlet fat extraction apparatus**—Complete with one flash with interchangeable ground glass joints made of Corning glass.

	Cap. of Flask	Cap. of Extractor	Price
(a)	100 ml	60 ml	Rs. 1,210
(b)	250 ml	100 ml	Rs. 1,390
(c)	500 ml	200 ml	Rs. 2,060
(d)	1000 ml	400 ml	Rs. 2,820
(e)	2000 ml	600 ml	Rs. 4,380

DA-4 **Tissue Homogenizers**—Hand operated with teflon or glass pestle

	2ml	5ml	10ml	15ml	25ml
Rs.	300	340	380	400	450

DA-5 **Craft's Electrically heated Distillation Apparatus** with continuous water flowing arrangement complete with metal stand and electric board. All glass borosilicate


Capacity	Single stage	Double stage	Triple stage
3 liter/hour	11,500	30,000	42,000
5 liter/hour	14,500	31,200	46,800

DA-8 **Craft Kjeldhal Distillation assembly**—complete with R.B. flask 1000 ml, socket, adopter, top funnel, splash head, condenser etc. **Rs. 2,400**


DA-9 **Craft's LANDSBURGER** Boiling apparatus **Rs. 1,450**

DA-10 **Specific Gravity Bottle** capacity

	5 ml	10 ml	25 ml	50 ml
Price each	Rs. 162	170	178	195




DA-1



DA-3


DA-4




DA-2




DA-10



DA-9




DA-5



DA-8


Biocraft's Glass Apparatus Made of Borosilicate Glass

DA-11	Weighing Bottle capacity				
	Capacity	5 ml	15 ml	25 ml	60 ml
	Price each	Rs. 175	188	212	275
DA-12	Pyknometer Double bend capacity				
	Capacity	5 ml	10 ml	20 ml	25 ml
	Price each	Rs. 165	181	193	210
DA-13	Wash Bottle with interchangeable Stopper				
	Capacity	250 ml	500 ml		
	Price each	Rs. 500	650		
DA-14	Gas washing bottle				
	Capacity	125 ml	250 ml	500 ml	
	Price each	Rs. 750	850	910	
DA-15	Craft's Condensers				
	Capacity	200 ml	300 ml	400 ml	
	Air Condenser	—	—	240	
	Leibig	412	475	525	
	Allihin	506	555	660	
	Coil	562	640	—	
DA-16	Chromatography Apparatus —It contains Buchner flask 250 ml, Tapadopter, column 200×18 mm with sintered disc, Reduction adopter, Reservoir 250 ml				Rs. 3,100
DA-17	Sprayer for Chromatography				
	Capacity	25 ml	50 ml	100 ml	
	Price each	Rs. 525	545	595	
DA-18	Nessler's Cylinder		50 ml	100 ml	
	Price each	Rs. 150	212		
DA-20	Kjeldhal Flask -Long Neck				
	Capacity	100 ml	300 ml	500 ml	
	Price each	Rs. 275	460	688	
DA-30	Gooch Crucible sintered				
	Capacity	15 ml	30 ml	50 ml	
	Price each	Rs. 140	160	225	
DA-31	Buchner Filtr with disc				
	Capacity	35 ml	80 ml	200 ml	500 ml
	Price each	Rs. 262	470	590	1,860
DA-32	Stalgnometer Stright—	Rs. 160,			
	Bent—	Rs. 180			
DA-33	Blood sugar tube capacity 25 ml.				Rs. 180 each
DA-34	Melting point tube 150×18 mm				Rs. 156 each
DA-35	Calcium chloride tube				Rs. 140 each
DA-36	Victor Meyer's apparatus				Rs. 1,350 each
DA-37	Viscometer College Pattern				Rs. 152 each



Biocraft's Glass Apparatus Made of Borosilicate Glass

DA-38	Biurettes with stop cock (Straight)					
	Cap.	Division	Borosilicate Glass			
	10 ml	1/20	Rs. 115			
	25 ml	1/10	Rs. 120			
	50 ml	1/10	Rs. 140			
	100 ml	1/5	Rs. 238			
DA-39	Biurettes with Teflon key (Rota flow)					
	Capacity	25 ml	50 ml	100 ml		
	Borosilicate Glass	Rs. 110	Rs. 140	Rs. 195		
DA-40	Biurettes with pinch cock					
	Cap.	Division	Borosilicate Glass			
	25 ml	1/10	Rs. 90			
	50 ml	1/10	Rs. 90			
	100 ml	1/5	Rs. 140			
DA-42	Pipettes : Graduated :					
	Capacity	Division	Borosilicate Glass			
	0.1 ml	0.01	Rs. 45			
	0.2 ml	0.01	Rs. 45			
	0.5 ml	0.02	Rs. 45			
	1 ml	0.1	Rs. 45			
	2 ml	0.1	Rs. 45			
	5 ml	0.1	Rs. 52			
	10 ml	0.1	Rs. 63			
	25 ml	0.2	Rs. 88			
DA-43	Pipettes : Single mark volumetric					
	Capacity	1 ml	2 ml	5 ml	10 ml	20 ml 25 ml
	Borosil Glass	Rs. 42	42	42	42	50 50
DA-44	Pasture Pipette —Hard Glass with Teat—	Rs. 30 each				
DA-45	Craft Filter Apparatus with Buchner Funnel					
	(A) Flask 250 ml, Funnel 80 ml	Rs. 990				
	(B) Flask 500 ml, Funnel 200 ml	Rs. 1,450				
DA-46	Crafts millipore type membrane filter holder assembly					
	Funnel Flask—	13 mm 250 ml flask	Rs. 3,650			
		25 mm 500 ml flask	Rs. 4,860			
DA-47	Milk Butyrometer —10% or 8% I.S.I.	Rs. 62				



Biocraft's Amber Colour Glass Ware

GAM-1	Reagent Bottle Narrow Mouth Flat Stopperd Dust Proof					
		30 ml	60 ml	125 ml	250 ml	500 ml 1000 ml
		Rs. 51	Rs. 58	Rs. 66	Rs. 99	Rs. 130 Rs. 230 each
GAM-2	Reagent Bottle Wide Mouth Flat Stopperd Dust Proof					
		125 ml	250 ml	500 ml		
		Rs. 83	Rs. 134	Rs. 195		
GAM-3	Dropping Bottles T.K. Pattern					
		30 ml	60 ml	125 ml		
		Rs. 97	Rs. 110	Rs. 144		