

Imported Cover Slips and Slides

Optically flat sheet thickness 0.13 to 0.15 mm.

These Cover Slips are very economical as each packet contains 100 Coverslip (Not 10 gm)—The weight of the 20×60 mm coverslip is approximately—58 gm. (nearly 6 times than the Inidan packing).

Imported square Coverslip No. 1

ICG-1	Size—18 × 18 mm	Per pkt. of 100 Nos. (Approx. 12 gm)	Rs. 56
ICG-2	Size—22 × 22 mm	Per pkt. of 100 Nos. (Approx. 20 gm)	Rs. 68
ICG-3	Size—18 mm	Per pkt. of 100 Nos. (Approx. 10 gm)	Rs. 77
ICG-4	Size—22 mm	Per pkt. of 100 Nos. (Approx. 20 gm)	Rs. 88

Imported rectangular Coverslip No. 1

ICG-5	Size—22 × 40 mm	Per pkt. of 100 Nos. (Approx. 36 gm)	Rs. 80
ICG-6	Size—22 × 50 mm	Per pkt. of 100 Nos. (Approx. 46 gm)	Rs. 84
ICG-7	Size—22 × 60 mm	Per pkt. of 100 Nos. (Approx. 58 gm)	Rs. 96
ICG-8	Size—20 × 30 mm	Per pkt. of 100 Nos. (Approx. 33 gm)	Rs. 80

Imported Glass Slides


Ground, polished edges thickness 1 mm. The slides are free from striation, scratch or dust as each slide is provide with tissue paper packing. Ready to use, most suitable for research, Pathology and Bacteriological studies. Size 25.4 × 76.2 (3"×1") thickness 1 mm-1.2 mm
Price **Rs. 62.00** per pkt. of 50 Nos.



Biocrafts Laboratory Dissecting Instruments

BDI-1 **Craft's dissecting trays with wax**—(Best quality). These special galvanized wax lined tray with locking walls with a deep layer of special tough wax which takes dissecting pins without cracking or chipping.

Size	Each Rs.	Size	Each Rs.
8"×6"×2.5"	195	14"×10"×3"	370
10"×8"×3"	280	14"×12"×3"	410
12"×10"×3"	298	15"×12"×3"	500




BDI-1

BDI-4 **Paraffin wax**—Milky white wax filling in dissection tray. Superior quality **Rs. 272** per Kg.

Biocraft'S Dissection Boxes & Instruments

BDI-9 **Craft's Dissection boxes**—(Best quality) with stainless steel instrument including razor
 Routine quality Special Delux quality
 (A) Set of 13 instruments **Rs. 185** **Rs. 380**
 (B) Set of 17 instruments **Rs. 250** **Rs. 530**
 (C) Set of 8 instruments (Plastic Pouch) **Rs. 180**



BDI-9

Biocrafts Laboratory Dissecting Instruments

(D) Medical grade Dissection Box-Superior Instruments **Rs. 780**

- BDI-11 **Scalpel Stainless Steel Superior**—1.5" blade **Rs. 28** each
 2" blade **Rs. 40** each
- BDI-11A **B.P. handle w/o blade dull polish 3" or 4"** **Rs. 30** each
- BDI-11B **Spare blade** **Rs. 64** Dozen
- BDI-12 **Forceps (Pointed or blunt) Stainless Steel**
 Size 4" 5" 6" 8" 10" 12" 15" 18"
Rs. 14 28 35 56 84 114 140 160
- BDI-13 **Scissors (Straight) Stainless steel**
 Size 4" 5" 6" 8"
Rs. 66 76 110 160
- BDI-14 **Scissors (Curved)—** 4" **Rs. 66** each
- BDI-15 **Bone Cutter Stainless Steel** 5" **Rs. 170** each
 6" **Rs. 185** each
- BDI-16 (a) Teasing Needles (all metal) **Rs. 16** each
 (b) Teasing Needles Plastic Handle **Rs. 5** each
- BDI-18 **Section cutting razor College Grade Superior** **Rs. 160** each
- BDI-21 **Spatula S. S. (with handle) 6" blade** **Rs. 80** each
 BDI-21 —do— 8" blade **Rs. 90** each




BDI-11



BDI-13

BDI-16



BDI-18



BDI-22



BDI-23

Craft'S Fine And Micro Dissection Forceps (delux) Stainless Steel

- BDI-22 Tissue forceps **Rs. 33**
- BDI-23 Fine tissue forceps **Rs. 33**
- BDI-24 Curved fine points forceps (Twiser) **Rs. 40**
- BDI-25 Fine points forceps **Rs. 36**
- BDI-26 Micro-dissection forceps **Rs. 36**
- BDI-27 Utility forcep curved **Rs. 36**



BDI-24




BDI-26



BDI-27

Craft'S Superior Delux Dissection Scissors (stainless Steel)

- BDI-28 Lister Curved Scissor 6" **Rs. 90**
- BDI-29 Mayo-Hedar Needle Holder 6" **Rs. 102**
- BDI-30 Operating Scissor 5" **Rs. 76**



BDI-28

Craft'S Common Surgical items

- BDIS-40 **Artery Forceps.** Stright or curved
 Size 5"—**Rs. 120** 6"—**Rs. 130** 8"—**Rs. 220** each
- BDIS-41 **Mosquito Artery Forcep.** Straight or curved
 Size 5"—**Rs. 120** 6"—**Rs. 130** each
- BDIS-42 **Needle Holder**
 Size 6"—**Rs. 140** 8"—**Rs. 230** each
- BDIS-43 **Allice Forceps**
 Size 6"—**Rs. 170** 8"—**Rs. 210** each
- BDIS-44 **Dissecting Forceps.** Plain or with teeth
 Size 5"—**Rs. 54** 6"—**Rs. 60** 8"—**Rs. 82** each
- BDIS-45 **Stainless Steel Tongs**
 Size 6"—**Rs. 85** 8"—**Rs. 100** 10"—**Rs. 130** 12"—**Rs. 170**



BDI-31

BDI-33

BDI-32